


No Party Like an ATAC-seq Party

Kevin Benac, Nima Hejazi, Hector Roux de Bézieux

Computational Biology 293
Spring 2018, Prof. Nir Yosef
21 February 2018


Introduction

- DNA presents different levels of compaction.
- Basic level : nucleosomes -- basic units of DNA packaging in eukaryotes cells.
- Double helix wrapped around 8 histone protein cores.



Introduction

- Then folded through a series of successively higher order structures to eventually form a chromosome.
- Both compacts DNA and creates an additional layer of regulatory control, which ensures correct gene expression.



ATAC-seq

- Recent 2-step technique, Buenrostro et al. (2013), that aims at identifying accessible regions of the chromatin + nucleosomal position.
- Relies on the action of mutated Tn5 transposase.
- Enzyme that catalyzes the movement of transposons to other parts in the genome.



ATAC-seq

Buenrostro et al. used these method on B-cell lines to


- identify regions of open chromatin;
- identify nucleosome-bound and nucleosome-free positions in regulatory regions;
- infer the positions of DNA binding proteins

Applied to assay personal T-cell epigenome of a healthy volunteer via standard serial blood draws

ATAC-seq

Proved to be more sensitive and more performant to measure and interpret the epigenome compared to FAIRE-seq and DNase-seq


- Faster;
- Requires fewer cells (500 to 50,000).



ATAC-seq probes chromatin accessibility with transposons

- Uses Tn5 transposase to integrate adapter into accessible chromatin regions
- Steric hindrance makes transposition of less accessible chromatin less probable
- Amplifiable DNA fragments suitable for high-throughput sequencing are generated at open chromatin locations
- Simple two-step process involving Tn5 transposase insertion and PCR for amplification

ATAC-seq probes chromatin accessibility with transposons



ATAC-seq probes chromatin accessibility with transposons

- Accurate, sensitive measure of chromatin accessibility
- ATAC-seq has similar signal-to-noise as DNase-seq, even when the latter is prepared with many more cells
- ROC curves display similar sensitivity and specificity for both ATAC-seq and DNase-seq
- Correlate well with active chromatin, not just regions of transposase preference
- Peak intensities highly correlated between ATAC-seq, DNase-seq

ATAC-seq insert sizes disclose nucleosome positions


- Paired-end reads provide a wealth of information about nucleosome packing and positioning
- High-resolution readout of nucleosome-associated and nucleosome-free regions in regulatory elements genome-wide

ATAC-seq insert sizes disclose nucleosome positions


- Different functional states of chromatin have differing accessibility “fingerprints” -- may be read with ATAC-seq
- Reveals differentially accessible forms of chromatin, hypothesized to exist *in vivo* but difficult to confirm

ATAC-seq reveals patterns of nucleosome - Transcriptome factor (TF) spacing

- Use CHIP-Seq to find where TF binds
- Compare nucleosome position to TF position
- Use unsupervised hierarchical clustering



ATAC-seq reveals patterns of nucleosome - Transcription factor (TF) spacing



ATAC-seq reveals patterns of nucleosome - Transcriptome factor (TF) spacing


“The interplay between precise nucleosome positioning and locations of DNA binding factor immediately suggests specific hypotheses for mechanistic studies, a potential advantage of ATAC-seq.”

ATAC-seq footprints infer factor occupancy genome-wide

Assumption: DNA-sequences occupied by DNA-binding proteins are protected for transposition

Results


One locus



ATAC-seq footprints infer factor occupancy genome-wide


Results

Average over all loci















Possible Personal Epigenomics?

ATAC-Seq enables epigenomic analysis on clinical timescales



Thank you! Questions?



- Body copy Lorem ipsum dolor sit amet, consectetur adipiscing elit, sed do eiusmod tempor incididunt ut labore et dolore magna aliqua.

- Body copy Lorem ipsum dolor sit amet, consectetur adipiscing elit, sed do eiusmod tempor incididunt ut labore et dolore magna aliqua.

- M

- Por

- incididunt ut labore et dolore magna aliqua.

- Body copy Lorem ipsum dolor sit amet, consectetur adipiscing elit, sed do eiusmod tempor incididunt ut labore et dolore magna aliqua.

Why are you reading this? This is clearly not fit for human consumption

Hope you enjoyed the talk!